


# INOGEN ONE<sup>G5™</sup> System



Inogen One  
AC Power Supply\*  
Item #BA-501



Inogen One G5  
DC Power Cable\*  
Item #BA-306



Inogen One G5  
Lithium Ion Battery\*  
Up to 6.5 hours run time  
Item #BA-500



Inogen One G5  
Carry Bag\*  
Item #CA-500



Inogen One G5  
Backpack\*\*  
Item #CA-550



Inogen One<sup>®</sup>G5  
Extended Life Lithium Ion Battery\*\*  
Up to 13 hours run time  
Item #BA-516

## The Inogen One G5 System (IS-500):

- Compact and Lightweight
- Flow Settings 1-6
- Intelligent Pulse Oxygen Delivery
- Designed for 24 Hour Use
- Easy-to-read LCD display
- Inogen Connect App

\*Included with system

\*\*Accessory sold separately

## Oxygen. Anytime. Anywhere®

Once again, Inogen innovation sets a new standard in oxygen delivery.

The Inogen One G5 offers the most oxygen per pound for a portable oxygen concentrator on the market today. With flow settings from 1-6, the Inogen One G5 is designed to dramatically increase independence for most users 24/7. The Inogen One G5's user interface has been designed with an easy-to-read LCD display and simple control functions. It can be charged at home or on the go - giving you the freedom of Oxygen. Anytime. Anywhere.

The Inogen One G5 is your one solution, your single solution, for oxygen at home or away.

<b>Weight</b>	4.7 lbs. (includes single battery)
<b>Size</b>	Length: 7.19" Width: 3.26" Height: 8.15" (includes single battery)
<b>Oxygen Flow</b>	Pulse dose delivery system Six flow settings 1, 2, 3, 4, 5, 6
<b>Power</b>	AC Power Supply 100-240V, 50-60Hz (auto sensing to allow worldwide use) DC Power Cable: for mobile use in car
<b>Battery</b>	Duration (single battery): Up to 6.5 hours Recharge up to 3 hours with AC or DC power Duration (double battery): Up to 13 hours Recharge up to 6 hours with AC or DC power
<b>Noise Level</b>	38 dBA*
<b>Warranty</b>	3 year or lifetime limited warranty**
<b>Operation</b>	Simple control functions and easy-to-read LCD display Inogen Connect App
<b>Use</b>	Designed for 24/7 use and conforms to all applicable FAA requirements for POC use on board an aircraft

\*At flow setting 2

\*\*See Inogen One® Limited Warranty Statement for complete warranty information

Manufactured by:



326 Bollay Drive  
Goleta, California 93117  
www.inogen.com



Conec Care BV  
Hondiusstraat 28L  
6827 DE Arnhem  
the Netherlands

+31854899400  
info@conec.care  
www.conec.care

